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Input: AC 120~240V

Output: DC 12V/1A

Rechargeable battery: 6V/4AH (Optional)

POWER CONSUMPTION

Aprox. 12mA

Aprox. 36mA with backlight.

Battery hour per charge: approx. 320 hours (backlight off)

KEYS FUNCTION

To turn off the balance, keep pressed this key until display shows "Off" for 2 seconds.

Press this key to turn on the balance.

Function 1. To select the desired weight unit.
Function 2. To exit from setup mode.

Function 1. To reset the weight to zero "0", but the display value has to be lesser than $\pm 2\%$ of maximum capacity.
Function 2. To move one space to the left or downward in setup mode.

Function 1. To subtract the container's weight.
Function 2. To move one space to the right or upward in setup mode.

Function 1. To view gross or net weight when the balance is on tare status.
All other keys will be disabled when gross weight is activated.
Function 2. To increase values upward in setup mode.

Function 1. Key of confirmation in setup mode.

STANDARD CALIBRATION

Press the  and  key at the same time to access to standard calibration.

ECF-1	Weight calibration (Zero and Span)
ECF-2	Zero calibration
ECF-3	Span calibration

Weight calibration ECF-1

1. Press the  key to access to mode ECF-1 or press the  key to exit and return to normal weighing mode.
2. The display will show CALZ, press the  key.
3. Use the keys ◀, ▶, ▲, 0 ~ 9 to introduce the calibration weight value and confirm pressing the  key.
4. Place the calibration weight on the platter and press the  key.
5. The balance will automatically calibrate and return to normal weighing mode.

Zero calibration ECF-2

1. Press the  key to access to mode ECF-2 or press the  key to exit and return to normal weighing mode.
2. The display will show CALZ, press the  key.
3. The balance will automatically perform zero calibration and return to normal weighing mode.

Span calibration ECF-3

1. Press the  key to access to mode ECF-3 or press the  key to exit and return to normal weighing mode.
2. Use the keys ◀, ▶, ▲, 0 ~ 9 to introduce the calibration weight value and confirm pressing the  key.
4. Place the calibration weight on the platter and press the  key.
5. The balance will automatically perform span calibration and return to normal weighing mode.

FUNCTIONS SETTING

To access to functions setting, press the  and  key.

UF-1 Internal value / display voltage

1. Press the  key to view the internal value.
2. Press again the  key to view the battery voltage.
3. To exit and return to normal weighing press the  key.

UF-2 Weight limits

1. Press the  key to access to weight limits configuration.
2. The display will show the message "0000.0L"
3. Use the keys  and  to introduce the low limit value, with key  to increase the digit value.
4. Press the  key to confirm.
5. The display will show "0000.0H"
6. Use the keys  and  to introduce the high limit value, with key  to increase the digit value.
7. Press the  key to confirm.

UF-3 Auto-turn off

Modes:

- **AoFF 00** – Auto-turn off disable.
- **AoFF 01** – The balance will automatically turn off after 1 minute of non use.
-

The auto-turn off can be configured up to 99 minutes AoFF01~AoFF99, the default value is AoFF10 (10 minutes)

1. Press the  key to access to auto-turn off configuration.
2. Use the keys  and  to move from one digit to other and the key  to select the desired mode.
3. Press the  key to confirm.

UF-4 Backlight

Modes:

- **Lit Off:** Disable.
- **Lit A:** Automatic
- **Lit on:** Always activated.

Default value: Lit A (automatic)

1. Press the  to access to backlight configuration.
2. Use the  key to select the desired mode.
3. Press the  key to confirm.

UF-5 HOLD function (keeps the reading fixed on display for few seconds after removing the load from the platter)

Modes:

- **HOLD 0** : Disable.
- **HOLD 1** : Weighing animals mode.
- **PctXXX**: To configurate the range from 001~100 in HOLD 1 mode.
- **timEXX**: To configurate the time from 1 2 4 8 16 32 64 in HOLD 1 mode

Default value: HOLD 0

1. Press the  key to access to HOLD function configuration.
2. Use the  key to select the desired mode.
3. Press the  key to confirm.

UF-8 Zero configuration

- ZP 0 Disable function.
- ZP 1 1 division not displayed
- ZP 2 2 divisions not displayed
- ZP 3 3 divisions not displayed
- ZP 4 4 divisions not displayed
- ZP 5 5 divisions not displayed

- This function is locked if HOLD 1 function is activated.

Default value: ZP 0

1. Press the  key to access to configuration.
2. Use the  key to select the desired mode.
3. Press the  key to confirm.

UF-9 Gravity Value setup

Although the balance allows the Standard Calibration, the G value setup provides a greater accuracy. This pre-calibration is performed with the switch CAL in ADJ position. The value will be set to -00- and it will be replaced each time you enter a new value.

Gravity value of Distributor: configured before the weight calibration.
Gravity value of User: configured after the weight calibration.

The value of gravity will be denied if is greater than 9.83217 (Polar Gravity) or less than 9.78031 (Ecuador Gravity)

Default value: 9.79423

1. Press the  key to access to configuration.
2. The display will show the last value for 1 second.
3. Use the keys  and  to move from a digit to other and the key  to increase the digit value.
4. Press the  key to confirm

GUARANTEE

This scale is guaranteed for one year from the delivery date. The guarantee covers any fabrication defect of the material.

During this period GRAM PRECISION, covers the manpower and the spare parts necessary for the reparation of the scale.

This guarantee does not cover the failures caused by an inappropriate use or overcharge.

The guarantee does not cover the freight cost (transport) necessary to repair the scale.

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