



► **Large-surface plate**
Large-surface weighing plate, protected from lateral impact.



For your lab and quality control


- Lab scales and quality control.
- Incorporates high-resolution aluminium sensor.
- Backlit LCD display with 21-mm digits.
- Underside weighing for density calculation.
- Large-size weighing plate.
- Keypad with basic functions.

► **Functions:**

- Cumulative tare of receptacles and products.
- Calculation of number of pieces.

External dimensions (mm)

EH



Model	EH-500	EH-1000	EH-2000	EH-3000	EH-6000
Reference	70542	70543	70544	70545	70546
Capacity	500 g	1,000 g	2,000 g	3,000 g	6,000 g
Resolution	0.01 g	0.01 g	0.01 g	0.1 g	0.1 g
Linearity	+/- 0.02 g	+/- 0.025 g	+/- 0.03 g	+/- 0.15 g	+/- 0.2 g
Reproducibility	+/- 0.01 g	+/- 0.015 g	+/- 0.02 g	+/- 0.1 g	+/- 0.15 g
Eccentricity	+/- 0.02 g	+/- 0.025 g	+/- 0.03 g	+/- 0.25 g	+/- 0.3 g
Reading time	less than 1.5 sec.				
Power	220V / 50-60 Hz				
Working temperature	+5°C / +40°C				
Warm-up time	10 minutes				
Calibration	external				
Weighing units	grams, ounces, pieces				
Air humidity	max 80% without condensation				
Hydrostatic weighing	underside bracket for hook				
Dim. of weighing plate (mm)	133 x 182		143 x 192		
Product dim. (mm)	225 x 220 x 75		225 x 220 x 80		
Product net weight (g)	1.2				
Packaging dim. (mm)	305 x 245 x 110				
Total weight with packaging (kg)	1.5				
RRP (€)	---	---	---	---	---