



► **Large-surface plate**
Large-surface weighing plate, protected from lateral impact.



Large weighings with precision


- High-resolution, high-capacity scales.
- Incorporates high-resolution aluminium sensor.
- Underside weighing for density calculation.
- LED display with 22-mm digits.
- Large-size weighing plate.
- RS232C data output.

► **Functions:**

- Cumulative tare of receptacles and products.
- Calculation of number of pieces.

External dimensions (mm)

EM



| Model | EM-10K | EM-20K | EM-30K |
|----------------------------------|------------------------------|------------|-----------|
| Reference | 70876 | 70877 | 70878 |
| Capacity | 10 kg | 20 kg | 30 kg |
| Resolution | 0.1 g | 0.1 g | 0.1 g |
| Linearity | +/- 0.2 g | +/- 0.2 g | +/- 0.2 g |
| Reproducibility | +/- 0.1 g | +/- 0.15 g | +/- 0.2 g |
| Eccentricity | +/- 0.2 g | +/- 0.25 g | +/- 0.3 g |
| Reading time | less than 1.5 sec. | | |
| Power | 220V / 50-60 Hz | | |
| Working temperature | +5°C / +40°C | | |
| Warm-up time | 10 minutes | | |
| Calibration | external | | |
| Weighing units | grams, ounces, pieces | | |
| Interface | RS232C | | |
| Air humidity | max 80% without condensation | | |
| Hydrostatic weighing | underside bracket for hook | | |
| Dim. of weighing plate (mm) | 340 x 240 | | |
| Product dim. (mm) | 375 x 390 x 130 | | |
| Product net weight (g) | 4.5 kg | | |
| Packaging dim. (mm) | 430 x 455 x 190 | | |
| Total weight with packaging (kg) | 5.5 kg | | |
| RRP (€) | --- | --- | --- |