

# GRAM WHITE MEL



## Large surface weighing pan

Large surface weighing pan with low profile.




## Extra-flat baby scales

- Low platform height of only 40 mm.
- Large weighing tray, the entire surface is usable.
- Tare function for baby clothes.
- Very resistant ABS structure.
- Automatic self-disconnection.
- Display with illumination for a perfect vision.
- Capacity of 20 kg and 5 grams of resolution.

### External dimensions (mm)

WHITE MEL



Model	WHITE MEL
Reference	4451
Capacity	20 kg
Resolution	10 g
Alimentación	100-240V_50/60Hz
Power	4 Batteries AA 1,5V.
Working temperature	0°C /+40°C
Calibration	external
Weighing units	g (kg), lb, oz
Dim. of weighing plate (mm)	550 x 260
Product dim. (mm)	550 x 330 x 40
Product net weight (kg)	1,6 kg
Packaging dim. (mm)	555 x 335 x 50
Total weight with packaging (kg)	2,1 kg